








Anterior Lumbar Articulating Arm System



ANTERIOR LUMBAR ARTICULATING ARM SYSTEM #SL82048

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
Tray 1					
A		43903AC	Elite III Rail Clamp with 2 Cam Joints 16"	2	Attaches to the rail of the OR table and holds the Articulating Arms in place.
B		SL42191G	Quick Angle Articulating Arm 12" x 21"	2	Inserts into Rail Clamp and holds retractor blades in place.
		60020	T-Handle	1	May be used on Articulating Arm to angle retractor blades 45 degrees up or down.
		42228	Manipulator Handle	1	May be used with the Articulating Arm to make controlled adjustments after blade has been placed.
		50000G	Instrument Case (23" x 11" x 3 1/2")	1	Holds instruments for sterilization, storage, and transportation.

ANTERIOR LUMBAR ARTICULATING ARM SYSTEM #SL82048 (Continued)

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
Tray 2					
		43902ACL	Elite III Rail Clamp with 1 Cam Joint 16"	1	Attaches to the rail of the OR table and holds the Articulating Arms in place.
		SL42191G	Quick Angle Articulating Arm 12" x 21"	2	Inserts into Rail Clamp and holds retractor blades in place.
		60020	T-Handle	1	May be used on Articulating Arm to angle retractor blades 45 degrees up or down.
		42228	Manipulator Handle	1	May be used on Articulating Arm to angle retractor blades 45 degrees up or down.
		SL46432ET	Radio Richardson 50mm x 160mm (2" x 6 1/4")	1	Abdominal wall retraction for procedures requiring a larger incision.
		SL46434ET	Radio Richardson 50mm x 180mm (2" x 7")	1	
		SL46436ET	Radio Richardson 50mm x 200mm (2" x 7 7/8")	1	
		50000G	Instrument Case (23" x 11" x 3 1/2")	1	Holds instruments.
Tray 3					
		SL46558	Radio Concave Short Neck Renal 25mm x 80mm (1" x 3 1/8")	4	Renal vein or vessel retraction.
		SL46560	Radio Concave Short Neck Renal 25mm x 100mm (1" x 4")	4	
		SL46562	Radio Concave Short Neck Renal 25mm x 120mm (1" x 4 3/4")	4	
		SL46564	Radio Concave Short Neck Renal 25mm x 140mm (1" x 5 1/2")	4	
		SL46566	Radio Concave Short Neck Renal 25mm x 160mm (1" x 6 1/4")	4	
		SL46568	Radio Concave Short Neck Renal 25mm x 180mm (1" x 7")	4	
		SL46570	Radio Concave Short Neck Renal 25mm x 200mm (1" x 7 7/8")	4	
		50000L	Instrument Tray - Lid	1	Tray Lid
		50000RV	Instrument Tray - Renal Vein Blades	1	Drops into Tray Base. Organizes blades.
		50000U1	Instrument Tray - Base	1	Tray Base