








Bariatric Liver Transplant



BARIATRIC SYSTEM #SL82001

Non S-Lock system also available, may require a longer lead time.

| REF | IMAGE | PART # | ITEM DESCRIPTION | QTY | APPLICATION |
|---------------|---|-------------|---|-----|---|
| Tray 1 | | | | | |
| A |  | 43905AC | Elite III Rail Clamp w/2 Cam Joints 22" | 2 | Attaches to the rail of the OR table and holds the bilateral crossbar and angled arms in position. |
| B |  | 41900PHL | 41 1/2" Bilateral Crossbar Hinged (14" x 10" x 14") | 1 | Attaches to the rail clamps to provide bilateral stability, and holds the retractor handles in place around the incision. |
| C |  | 44124W | 24" Angled Arm (11" x 13") @ 45° | 2 | Attaches to the crossbar and holds the retractor handles in place around the incision. |
| | | 50000BL | Instrument Case (26" x 10" x 5") | 1 | Holds instruments. |
| Tray 2 | | | | | |
| D |  | SL42126WGP | Cam II Clip-on Quick Angle 8" | 5 | Attaches retractor blades to frame and angles retractor blades as needed, 45° up or down. |
| |  | SL45006CGP | Micro-Adjustable II Clip-on Quick Angle 10" | 2 | Attaches blades to frame and angles blades as needed, 45° up or down. |
| E |  | SL45006CGRP | Micro-Adjustable II Clip-on Quick Angle 15" | 2 | Provides precise micro-adjustments. |

BARIATRIC SYSTEM #SL82001 (Continued)

| REF | IMAGE | PART # | ITEM DESCRIPTION | QTY | APPLICATION |
|-----|---|----------|--|-----|---|
| |  | 60020 | T-Handle | 2 | May be used on Quick Angle Handle to angle blade 45° up or down. |
| |  | SL46119C | Concave Round Neck Renal 25mm x 190mm (1" x 7 1/2") | 2 | Renal vein or vessel retraction. |
| |  | SL46122A | Malleable 51mm x 254mm (2" x 10") | 1 | Soft tissue retraction. |
| | | SO46124A | Malleable 76mm x 254mm (3" x 10") | 1 | |
| | | SO46124B | Malleable 102mm x 254mm (4" x 10") | 1 | |
| |  | SO46130 | Kelly 76mm x 89mm (3" x 3 1/2") | 2 | Abdominal or incisional wall retraction. |
| |  | SO46136A | Fence 102mm x 203mm (4" x 8") | 1 | Lateral retraction of abdominal wall on large procedures. |
| |  | SO46141 | Balfour 83mm x 73mm (3 1/4" x 3") | 2 | Costal margin or abdominal wall retraction. |
| | | SO46141D | Balfour 83mm x 123mm (3 1/4" x 5") | 2 | |
| | | SO46141F | Balfour 83mm x 175mm (3 1/4" x 6 7/8") | 1 | |
| |  | SO46153 | Richardson 51mm x 178mm (2" x 7") | 1 | Abdominal wall retraction. |
| |  | SL46161 | Harrington 114mm x 254mm (4 1/2" x 10") | 1 | Organ or soft tissue retraction. |
| |  | SO46170C | Non-Slip Balfour R 77mm x 94mm (3" x 3 3/4") | 1 | Costal margin retraction. Aligns with anatomy, reducing slippage. |
| | | SO46171C | Non-Slip Balfour L 77mm x 94mm (3" x 3 3/4") | 1 | |
| |  | SL46191 | Malleable Finger 152mm (6") | 1 | Gentle organ retraction. |
| |  | SL46146 | Splanchnic 51mm x 152mm (2" x 6") | 1 | Soft tissue retraction. |
| | | 50000G | Instrument Case (23" x 11" x 3 1/2") | 1 | Holds instruments. |