



## Sternotomy \*













\* Shown with Morse Sternal Spreader ( #300251 )

## Thoracotomy \*\*



\*\* Shown with thoracotomy spreader, not sold by Thompson Surgical.

**BOLLING MITRAL VALVE RETRACTOR SYSTEM #91072CV**

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
<b>Tray 1</b>					
		43100	Malleable Paddle 5" (127mm)	1	Malleable atrial or septal retraction
		43101	Malleable Paddle 6 1/2"	1	
		43102	Malleable Paddle 8" (203mm)	1	
		43103	Cushing 5" (127mm)	1	Atrial retraction during transseptal approach.
		43104	Cushing 6 1/2" (165mm)	1	
		43105	Cushing 8" (203mm)	1	
		43106	Small Basket 5" (127mm)	1	Left atrial wall retraction.
		43107	Small Basket 6 1/2" (165mm)	1	
		43108	Small Basket 8" (203mm)	1	
		43109	Medium Basket 5" (127mm)	1	Left atrial wall retraction.
		43110	Medium Basket 6 1/2" (165mm)	1	
		43111	Medium Basket 8" (203mm)	1	
		43112	Large Basket 5" (127mm)	1	Left atrial wall retraction.
		43113	Large Basket 6 1/2" (165mm)	1	
		43114	Large Basket 8" (203mm)	1	
A		47204	Universal Joint	3	Attaches retractors to CV Clamps.
		47202N	CV Clamp L-Bar	1	Mounts to sternal or rib retractors and holds multiple Universal Joints.
B		47200N	CV Clamp Mini L-Bar	2	Mounts to sternal or rib retractors and holds one Universal Joint.
C		47201N	CV Clamp T-Bar	1	Mounts to sternal or rib retractors and holds multiple Universal Joints.
		50000BOL	Instrument Case - Bolling Mitral Valve (23" x 10 3/4" x 2 1/2")	1	Organizes and holds instruments for sterilization, storage, and transportation.