


OneFrame™ Frame Kit



ONEFRAME™ FRAME KIT #SL82063

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
Tray 1					
A		43902AC	Elite III Rail Clamp with 1 Cam Joint 18"	1	Attaches to the rail of the OR table and holds the crossbar in position.
B		47500	OneFrame®	1	Attaches to the rail clamp and holds the retractor handles in place around the incision.
		44013H	Extension Arm 18" with Spring Hinge	1	Attaches to distal end of Angled Armsto provide extra reach, closure, or caudal retraction.
		44118B:	Angled Closure Arm 18" (9" x 9")	2	Attaches to distal ends of frame arms to create enclosed frame.
C		SL42126WGP	Cam II Clip-On Quick Angle 8"	5	Attaches retractor blades to the OneFrame and angles retractor blades as needed, 45° up or down.
		SL45006CGP	Micro-Adjustable II Clip-On Quick Angle 10"	1	Attaches blades to frame and angles blades as needed, 45° up or down. Provides precise micro-adjustments.
		60020	T-Handle	2	May be used on Quick Angle Handle to angle blade 45° up or down.
		50000G	Instrument Case (23" x 11" x 3 1/2")	1	Holds instruments.

Anterior Lumbar Blade Kit



ANTERIOR LUMBAR BLADE KIT #SL82065

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
Tray 1					
		SL46192ET	Radio. Reverse Lip Malleable 25mm x 200mm	2	Vessel or soft tissue retraction.
		SL46193ET	Radio. Reverse Lip Malleable 32mm x 200mm	2	
		SL46260RET	Radiolucent Reverse Lip 25mm x 100mm	2	Reverse Lip engages at the lateral aspect of the disc and vertebral body to expose the anterior lumbar spine.
		SL46261RET	Radiolucent Reverse Lip 25mm x 110mm	2	
		SL46262RET	Radiolucent Reverse Lip 25mm x 120mm	2	
		SL46263RET	Radiolucent Reverse Lip 25mm x 130mm	2	
		SL46265RET	Radiolucent Reverse Lip 25mm x 150mm	2	
		SL46267RET	Radiolucent Reverse Lip 25mm x 170mm	2	
		SL46269RET	Radiolucent Reverse Lip 25mm x 190mm	2	
		SL46271RET	Radiolucent Reverse Lip 25mm x 210mm	2	
		SL46273RET	Radiolucent Reverse Lip 25mm x 230mm	2	
		SL46275RET	Radiolucent Reverse Lip 25mm x 250mm	2	
		SL46280RET	Radiolucent Reverse Lip 32mm x 100mm	2	
		SL46282RET	Radiolucent Reverse Lip 32mm x 120mm	2	
		SL46284RET	Radiolucent Reverse Lip 32mm x 140mm	2	Reverse Lip engages at the lateral aspect of the disc and vertebral body to expose the anterior lumbar spine.
		SL46286RET	Radiolucent Reverse Lip 32mm x 160mm	2	
		SL46288RET	Radiolucent Reverse Lip 32mm x 180mm	2	
		SL46290RET	Radiolucent Reverse Lip 32mm x 200mm	2	
		SL46292RET	Radiolucent Reverse Lip 32mm x 220mm	2	
		SL46293RET	Radiolucent Reverse Lip 32mm x 230mm	2	
		SL46294RET	Radiolucent Reverse Lip 32mm x 240mm	2	Tapered design is ideal for working around osteophytes.
		SL46295RET	Radiolucent Reverse Lip 32mm x 250mm	2	
		SL46260TET	Radio. Reverse Lip Tapered 25mm x 100mm	2	
		SL46265TET	Radio. Reverse Lip Tapered 25mm x 150mm	2	
		SL46270TET	Radio. Reverse Lip Tapered 25mm x 200mm	2	
		SL46119BET	Radiolucent Malleable Concave Round Neck Renal 25mm x 140mm	2	Renal vein or vessel retraction.
		SL46119CET	Radiolucent Malleable Concave Round Neck Renal 25mm x 190mm	2	
	(see image below)	SL46564	Radio. Conc. Short Neck Renal 25mm x 140mm	4	Rigid, low profile vessel retraction.
		SL46569	Radio. Conc. Short Neck Renal 25mm x 190mm	4	
		SL46575	Radio. Conc. Short Neck Renal 25mm x 250mm	4	
		51234	Suction for Anterior Lumbar Surgery	1	Suction for disc space, as needed.