

Bilateral Anterior Cervical



CERVICAL SYSTEM WITH LITEPATH #SL82043

Non S-Lock and Titanium systems also available, may require longer lead times.

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
Tray 1					
A		43902ACS	Elite III Rail Clamp w/1 Cam Joint 11"	2	Attaches to the rail of the OR table and holds the angled arm in position.
B		44120	20" Angled Arm (10" x 10") @ 45°	2	Attaches to the rail clamp and holds the retractor handles in place around the incision.
C		SL42126WGP	Cam II Clip-on Quick Angle 8"	4	Attaches retractor blades to frame and angles retractor blades as needed, 45° up or down.
		SL45011CGP	Mini Micro-Adjustable II Clip-on Quick Angle 9"	1	Attaches blades to frame and angles blades as needed, 45° up or down. Provides precise micro-adjustments.
		60020	T-Handle	2	May be used on Quick Angle Handle to angle blade 45° up or down.
		42228	Manipulator Handle	2	May be used with a Thompson retractor handle to make controlled adjustments after blade has been placed and attached to frame.
		44012H	Extension Arm 12" with Snap Spring Joint	2	May be added to the angled arm if superior retraction is needed.
		50000G	Instrument Case 22" x 11" x 3 1/2"	1	Holds instruments.
		5844	Rail Extender 15" Single Clamp	2	Attaches to OR table rail to increase rail length or width and provide more attachment options for the rail clamps.

CERVICAL SYSTEM WITH LITEPATH #SL82043 (Continued)

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
Tray 2					
		SL46312EBL	Radiolucent Concave with LitePath 23mm x 35mm	2	Medial and lateral retraction for multi and single level procedures.
		SL46314EBL	Radiolucent Concave with LitePath 23mm x 40mm	2	
		SL46316EBL	Radiolucent Concave with LitePath 23mm x 45mm	2	
		SL46318EBL	Radiolucent Concave with LitePath 23mm x 50mm	2	
		SL46320EBL	Radiolucent Concave with LitePath 23mm x 55mm	2	
		SL46322EBL	Radiolucent Concave with LitePath 23mm x 60mm	2	
		SL46324EBL	Radiolucent Concave with LitePath 23mm x 65mm	2	
		SL46326EBL	Radiolucent Concave with LitePath 23mm x 70mm	2	
		SL46328EBL	Radiolucent Concave with LitePath 23mm x 75mm	2	
		SL46330EBL	Radiolucent Concave with LitePath 23mm x 80mm	2	
		SL45160EBL	Radio. Conc. Longitudinal with LitePath 23mm x 30mm	2	Longitudinal, medial, or lateral retraction.
		SL45161EBL	Radio. Conc. Longitudinal with LitePath 23mm x 40mm	2	
		SL45162EBL	Radio. Conc. Longitudinal with LitePath 23mm x 50mm	2	
		SL45163EBL	Radio. Conc. Longitudinal with LitePath 23mm x 60mm	2	
		SL45164EBL	Radio. Conc. Longitudinal with LitePath 23mm x 70mm	2	
		SL46510EBL	Radio. Concave Micro with LitePath 14mm x 35mm	2	Medial or lateral retraction. Ideal for small exposures and longitudinal retraction in narrow exposures.
		SL46512EBL	Radio. Concave Micro with LitePath 14mm x 40mm	2	
		SL46514EBL	Radio. Concave Micro with LitePath 14mm x 45mm	2	
		SL46516EBL	Radio. Concave Micro with LitePath 14mm x 50mm	2	
		SL46518EBL	Radio. Concave Micro with LitePath 14mm x 55mm	2	
		SL46520EBL	Radio. Concave Micro with LitePath 14mm x 60mm	2	Medial or lateral retraction.
		SL45643EBL	Radio. Concave Tapered with LitePath 23mm x 35mm	2	
		SL45644EBL	Radio. Concave Tapered with LitePath 23mm x 40mm	2	
		SL45645EBL	Radio. Concave Tapered with LitePath 23mm x 45mm	2	
		SL45646EBL	Radio. Concave Tapered with LitePath 23mm x 50mm	2	
		SL45647EBL	Radio. Concave Tapered with LitePath 23mm x 55mm	2	
		SL45648EBL	Radio. Concave Tapered with LitePath 23mm x 60mm	2	
		SL45649EBL	Radio. Concave Tapered with LitePath 23mm x 65mm	2	Retraction under the mandible on high level anterior cervical procedures.
		SL45650EBL	Radio. Concave Tapered with LitePath 23mm x 70mm	2	
		SL45750EBL	Radio. Midline Right with LitePath 70mm x 83mm	1	Retraction under the mandible on high level anterior cervical procedures.
		SL45751EBL	Radio. Midline Left with LitePath 70mm x 83mm	1	
		SL45666EBL	Radiolucent Odontoid with LitePath 20mm x 90mm	1	Anterior odontoid and longitudinal cephalad retraction.
		46519	Depth Gauge - Anterior Cervical [Anodized]	1	Measures depth of incision for blade selection.
		50000EAC	Instrument Case - Expanded Anterior Cervical	1	Drops into Tray Base and organizes blades
		50000U1	Instrument Tray - Base	1	Tray Base
		50000L	Instrument Tray - Lid	1	Tray Lid

OPTIONAL ADD-ON LIGHTING #TRA00001

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
		40040	LitePath Integrated Double 13mm	5	Inserts into blades for illumination of wound.
		ROVER-SLS-S-5	Rover Disposable Light Source	1	Battery-powered light source for LitePath. (5 per box)