


## Muscle Retractor ( Posterior )



### MUSCLE RETRACTOR ADD-ON KIT #SL90060

NOTE: Hip Frame required for use.

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
<b>Tray 1</b>					
A		SL45006CGP	Micro-Adjustable II Clip-on Quick Angle 10"	1	Attaches blades to frame and angles blades as needed, 45° up or down. Provides precise micro-adjustments.
		SL45006CGRP	Micro-Adjustable II Clip-on Quick Angle 15"	1	
		60020	T-Handle	2	May be used on Quick Angle Handle to angle blade 45° up or down.
		SL46410T	Hibbs with Teeth 25mm x 75mm (1" x 3")	2	Muscle retraction.
		SL46414T	Hibbs with Teeth 25mm x 115mm (1" x 4 1/2")	2	
		SL46418T	Hibbs with Teeth 25mm x 165mm (1" x 6 1/2")	2	
		50000G	Instrument Case (23" x 11" x 3 1/2")	1	Holds instruments for sterilization, storage, and transportation.