



Radiolucent Liver / Oncology / HPB



RADIOLUCENT LIVER / ONCOLOGY SYSTEM #SL82008

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
Tray 1					
A		43903AC	Elite III Rail Clamp with 2 Cam Joints 18"	2	Attaches to OR table, holds crossbar and arms.
B		41900PH	34" Bilateral Crossbar Hinged 11 1/2" x 8 1/2" x 11 1/2"	1	Attaches to the rail clamps to provide bilateral stability, and holds the retractor handles in place around the incision.
		41900F	Bilateral Crossbar 7 1/4" x 13 1/2" x 7 1/4"	1	
		44112SP	20" Lower Abdominal Bar (14" x 6")	1	Attaches to the rail clamps or crossbar and holds the retractor handles in place around the incision.
		50000BL	Instrument Case	1	Holds instruments.
Tray 2					
C		44118N	18" Angled Arm (8" x 10") @ 45°	1	Attaches to the crossbar and holds the retractor handles in place around the incision.
		44124W	24" Angled Arm (11" x 13") @ 45°	2	
		44124N	24" Angled Arm (8" x 16") @ 45°	1	
		43112C	Cam III Joint Serrated 1/2" x 1/2"	2	Connects angled arms together to allow for expansion of frame.
D		SL42126WGP	Cam II Clip-on Quick Angle 8"	6	Attaches retractor blades to frame and angles retractor blades as needed, 45° up or down.
E		SL45006CGP	Micro-Adjustable II Clip-on Quick Angle 10"	2	Attaches blades to frame and angles blades as needed, 45° up or down. Provides precise micro-adjustments.
		SL45006CGRP	Micro-Adjustable II Clip-on Quick Angle 15"	2	

RADIOLUCENT LIVER / ONCOLOGY SYSTEM #SL82008 (Continued)

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
		60020	T-Handle	2	May be used on Quick Angle Handle to angle blade 45° up or down.
		50000G	Instrument Case (23" x 11" x 3 1/2")	1	Holds instruments.
Tray 3					
		SL46115ET	Malleable 13mm x 64mm (1/2" x 2 1/2")	2	Soft tissue retraction.
		SL46116ET	Malleable 13mm x 152mm (1/2" x 6")	2	
		SL46117ET	Malleable 13mm x 203mm (1/2" x 8")	2	
		SL46118ET	Malleable 25mm x 203mm (1" x 8")	1	
		SL46118PET	Malleable 25mm x 104mm (1" x 4")	1	
		SL46122ET	Malleable 51mm x 203mm (2" x 8")	1	
		SL46124ET	Malleable 64mm x 254mm (2 1/2" x 10")	2	
		SO46124AET	Malleable 76mm x 254mm (3" x 10")	1	
		SO46124BET	Malleable 102mm x 254mm (4" x 10")	1	
			SO46129ET	Kelly 64mm x 76mm (2 1/2" x 3")	
	SO46130ET		Kelly 76mm x 89mm (3" x 3 1/2")	2	
		SO46140CEB	Balfour [small lips] 65mm x 72mm (2 1/2" x 2 3/4")	1	Costal margin or abdominal wall retraction.
		SO46141EB	Balfour [lips] 83mm x 70mm (3 1/4" x 2 3/4")	2	
		SO46141CEB	Balfour [small lips] 83mm x 70mm (3 1/4" x 2 3/4")	2	
		SO46141DEB	Balfour [lips] 83mm x 127mm (3 1/4" x 5")	2	
		SL46142EB	Balfour [lips] 32mm x 32mm (1 1/4" x 1 1/4")	4	
		SL46153ET	Richardson 51mm x 178mm (2" x 7")	1	Abdominal wall retraction.
		SL46160ET	Harrington 64mm x 152mm (2 1/2" x 6")	1	Organ or soft tissue retraction.
		SO46188ET	St. Marks [with lip] 64mm x 178mm (2 1/2" x 7")	1	Soft tissue retraction.
		SL46191ET	Malleable Finger 152mm (6")	2	Gentle organ retraction.
		SO46241AEB	Wide Balfour 83mm x 57mm (3 1/4" x 2 1/4")	1	Costal margin or abdominal wall retraction.
		SO46241BEB	Wide Balfour 100mm x 61mm (4" x 2 3/8")	1	
		SO46241CEB	Wide Balfour 114mm x 63mm (4 1/2" x 2 1/2")	1	
		SO46170CEB	Non-Slip Balfour L 77mm x 94mm (3" x 3 3/4")	1	Costal margin retraction. Aligns with anatomy, reducing slippage.
		SO46171CEB	Non-Slip Balfour R 77mm x 94mm (3" x 3 3/4")	1	
		50000G	Instrument Case (22" x 11" x 3 1/2")	1	Holds instruments.