


Colon-Rectal RingTrack® Retractor Set



RINGTRACK® FRAME KIT #SL91081

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
Tray 1					
A		47902ACL	RingTrack Elite III Rail Clamp with 1 Cam Joint 16 "	1	Attaches to the rail of the OR table and holds the crossbar in position.
B		47437	RingTrack Horizontal Bar 14 ¼ "	1	Attaches to the rail clamp and holds the ratcheted ring in place.
C		47412	RingTrack Oval Ring Medium 10 ½ x 14 ½ "	1	Ratcheted Ring that holds retractor handles, ratchets, and blades around the incision.
D		SL45005G	RingTrack Quick Angle Handle 8 "	6	Attaches to S-Lock retractor blades, and inserts into a ratchet to be attached to ring.
E		47440	RingTrack Ratchet	6	Attaches handles to frame.
		60020	T-Handle	2	May be used on Quick Angle Handle to angle blade 45° up or down.
		50000W	Instrument Case (24" x 16" x 3")	1	Holds instruments.

COLON-RECTAL BLADE KIT #SL91007L

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
Tray 1					
		SL46186	Hill-Ferguson 25mm x 64mm (1" x 2 1/2")	1	Exposure for hemorrhoid surgery.
		SL46187	Rectal Blade 13mm x 117mm (1/2" x 4 1/2")	2	Soft tissue retraction to expose bowel.
		SL46143	Fenestrated 44mm x 89mm (1 3/4" x 3 1/2")	2	Abdominal or incisional wall retraction.
D		SL46170	Deaver 25mm x 127mm (1" x 5")	2	Soft tissue retraction.
		SL46185	Hill-Ferguson 32mm x 76mm (1 1/4" x 3")	1	Exposure for hemorrhoid surgery.
E		SO46187A	St. Marks [no lip] 64mm x 178mm (2 1/2" x 7")	1	Soft tissue retraction.
		SO46188	St. Marks [with lip] 64mm x 178mm (2 1/2" x 7")	1	
		SO46189	St. Marks [with lip] 64mm x 140mm (2 1/2" x 5 1/2")	1	
		50000G	Instrument Case (23" x 11" x 3 1/2")	1	Holds instruments.