


Hysterectomy RingTrack® Retractor Set

Abdominal Hysterectomy RingTrack Set Up



Vaginal Hysterectomy RingTrack Set Up



RINGTRACK® FRAME KIT #SL91081

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
Tray 1					
A		47902ACL	RingTrack Elite III Rail Clamp with 1 Cam Joint 16"	1	Attaches to the rail of the OR table and holds the crossbar in position.
B		47437	RingTrack Horizontal Bar 14 1/4"	1	Attaches to the rail clamp and holds the ratcheted ring in place.
C		47412	RingTrack Oval Ring Medium 10 1/2 x 14 1/2"	1	Ratcheted Ring that holds retractor handles, ratchets, and blades around the incision.
D		SL45005G	RingTrack Quick Angle Handle 8"	6	Attaches to S-Lock retractor blades, and inserts into a ratchet to be attached to ring.
E		47440	RingTrack Ratchet	6	Attaches handles to frame.
		60020	T-Handle	2	May be used on Quick Angle Handle to angle blade 45° up or down.
		50000W	Instrument Case (24" x 16" x 3")	1	Holds instruments.

GYNECOLOGY BLADE KIT #SL91020GYN

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
Tray 1					
		SL46174	Deaver 51mm x 127mm (2" x 5")	1	Soft tissue retraction.
		SL46180	Heany 25mm x 102mm (1" x 4")	1	Abdominal or incisional wall retraction.
		SL46160	Harrington 64mm x 152mm (2 1/2" x 6")	1	Organ or soft tissue retraction.
		SL46663	Double Angled 13mm x 70mm (1/2" x 2 3/4")	2	Soft tissue retraction.
		SL46660	Double Angled 25mm x 76mm (1" x 3")	2	
		SL46118A	Malleable 25mm x 254mm (1" x 10")	2	Soft tissue retraction.
		SL46118P	Malleable 25mm x 102mm (1" x 4")	2	
		SL46118	Malleable 25mm x 203mm (1" x 8")	2	
		SL46122	Malleable 51mm x 203mm (2" x 8")	2	
		SO46135	Fence 102mm x 127mm (4" x 5")	1	Lateral retraction of abdominal wall on large procedures.
		SL46150	Richardson 51mm x 102mm (2" x 4")	1	Abdominal wall retraction.
		SL46143	Fenestrated 44mm x 89mm (1 3/4" x 3 1/2")	2	Abdominal or incisional wall retraction.
		50000G	Instrument Case (23" x 11" x 3 1/2")	1	Holds instruments.
		5844	Rail Extender 15" with 1 Power Rail Clamp	1	Attaches to OR table rail to increase rail length or width and provide more attachment options for the rail clamps.