



PRC Rigid Holder + Flexible Holder




PRC Articulating Arm + Flexible Holder



LAPAROSCOPIC SYSTEM #80030

REF	IMAGE	PART #	ITEM DESCRIPTION	APPLICATION
A		41950	Power Rail Clamp (PRC)	Attaches to the rail of the OR table and holds the laparoscopic holding arm(s) in position.
B		42133DN	Flexible Laparoscopic Holder for PRC	Unlimited range of motion and positioning with adjustable tension from flexible to rigid.
C		42139N	Rigid Laparoscopic Holder for PRC	Precise, rigid positioning for instruments requiring few adjustments or holding excess weight. Additional rigid holder options available.
D		42138L	Laparoscopic Articulating Arm for PRC	Holds firmly with a wide range of motion.
E		33181	Adjustable Scope Grip 5-6 / 10-11mm	Holds scope firmly while protecting shaft. Fits scopes 5-6 or 10-11mm.
F		42206	Nathanson Hook Small	Organ retraction.
		42207	Nathanson Hook Medium	
		42208	Nathanson Hook Large	

LAPAROSCOPIC SYSTEM #80030 (Continued)

REF	IMAGE	PART #	ITEM DESCRIPTION	APPLICATION
		42160N	Due Head 3/5/10mm with Knob	Holds Laparoscopic Instruments
		50000BL	Instrument Case (26" x 10" x 5")	Holds Instruments in System

Additional attachments available, not shown.