


OneFrame™ Frame Kit



ONEFRAME™ FRAME KIT #SL82063

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
Tray 1					
A		43902AC	Elite III Rail Clamp with 1 Cam Joint 18"	1	Attaches to the rail of the OR table and holds the crossbar in position.
B		47500	OneFrame®	1	Attaches to the rail clamp and holds the retractor handles in place around the incision.
		44013H	Extension Arm 18" with Spring Hinge	1	Attaches to distal end of Angled Armsto provide extra reach, closure, or caudal retraction.
		44118B:	Angled Closure Arm 18" (9" x 9")	2	Attaches to distal ends of frame arms to create enclosed frame.
C		SL42126WGP	Cam II Clip-On Quick Angle 8"	5	Attaches retractor blades to the OneFrame and angles retractor blades as needed, 45° up or down.
		SL45006CGP	Micro-Adjustable II Clip-On Quick Angle 10"	1	Attaches blades to frame and angles blades as needed, 45° up or down. Provides precise micro-adjustments.
		60020	T-Handle	2	May be used on Quick Angle Handle to angle blade 45° up or down.
		50000G	Instrument Case (23" x 11" x 3 1/2")	1	Holds instruments.

General / Vascular Blade Kit



GENERAL / VASCULAR BLADE KIT #SL82066

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
		SL46119B	Concave Round Neck Renal 25mm x 140mm (1" x 5 1/2")	2	Renal vein or vessel retraction.
		SL46119C	Concave Round Neck Renal 25mm x 190mm (1" x 7 1/2")	2	
		SL46122V	Malleable Slotted 51mm x 203mm (2" x 8")	1	Soft tissue retraction.
		SO46124AV	Malleable Slotted 76mm x 203mm (3" x 8")	1	
		SL46118	Malleable 25mm x 203mm (1" x 8")	1	Soft tissue retraction.
		SL46122	Malleable 51mm x 203mm (2" x 8")	2	
		SL46124	Malleable 64mm x 254mm (2 1/2" x 10")	1	
		SO46124A	Malleable 76mm x 254mm (3" x 10")	1	
		SL46125	Richardson 19mm x 25mm (3/4" x 1")	1	Muscle retraction for exposure with small or mini-incision.
		SO46129	Kelly 64mm x 76mm (2 1/2" x 3")	2	Abdominal or incisional wall retraction.
		SO46130	Kelly 76mm x 89mm (3" x 3 1/2")	1	
		SO46134	Fence 102mm x 102mm (4" x 4")	1	Lateral retraction of abdominal wall on large procedures.
		SO46135	Fence 102mm x 127mm (4" x 5")	1	
		SO46136	Fence 102mm x 152mm (4" x 6")	1	
		SO46140	Balfour [lips] 70mm x 73mm (2 3/4" x 3")	1	Costal margin or abdominal wall retraction.
		SO46140C	Balfour [sm lips] 66mm x 73mm (2 1/2" x 3")	2	
		SO46141	Balfour [lips] 83mm x 73mm (3 1/4" x 3")	1	
		SO46141C	Balfour [sm lips] 78mm x 73mm (3" x 3")	1	
		SL46143	Fenestrated 44mm x 89mm (1 3/4" x 3 1/2")	2	Abdominal or incisional wall retraction.
		SL46144	Fenestrated 48mm x 127mm (1 7/8" x 5")	2	

GENERAL / VASCULAR BLADE KIT #SL82066 (Continued)

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
		SL46145	Splanchnich 51mm x 114mm (2" x 4 1/2")	1	Soft tissue retraction.
		SL46146	Splanchnich 51mm x 152mm (2" x 6")	1	
		SL46151	Richardson 51mm x 127mm (2" x 5")	2	Abdominal wall retraction.
		SL46160	Harrington 64mm x 152mm (2 1/2" x 6")	1	Organ or soft tissue retraction.
		SL46174	Deaver 51mm x 127mm (2" x 5")	1	Soft tissue retraction.
		SO46176	Deaver 76mm x 127mm (3" x 5")	1	
		SO46188	St. Marks [lip] 64mm x 178mm (2 1/2" x 7")	1	Soft tissue retraction.
		SL46191	Malleable Finger 152mm (6")	1	Gentle organ retraction.
		SO46241A	Wide Balfour 83mm x 57mm (3 1/4" x 2 1/4")	1	Costal margin or abdominal wall retraction.
		SO46241B	Wide Balfour 100mm x 61mm (4" x 2 3/8")	1	
		SO46241C	Wide Balfour 114mm x 63mm (4 1/2" x 2 1/2")	1	
		50000G	Instrument Case (23" x 11" x 3 1/2")	1	Holds instruments.