


OneFrame™ Frame Kit



ONEFRAME™ FRAME KIT #SL82063

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
Tray 1					
A		43902AC	Elite III Rail Clamp with 1 Cam Joint 18"	1	Attaches to the rail of the OR table and holds the crossbar in position.
B		47500	OneFrame®	1	Attaches to the rail clamp and holds the retractor handles in place around the incision.
		44013H	Extension Arm 18" with Spring Hinge	1	Attaches to distal end of Angled Armsto provide extra reach, closure, or caudal retraction.
		44118B:	Angled Closure Arm 18" (9" x 9")	2	Attaches to distal ends of frame arms to create enclosed frame.
C		SL42126WGP	Cam II Clip-On Quick Angle 8"	5	Attaches retractor blades to the OneFrame and angles retractor blades as needed, 45° up or down.
		SL45006CGP	Micro-Adjustable II Clip-On Quick Angle 10"	1	Attaches blades to frame and angles blades as needed, 45° up or down. Provides precise micro-adjustments.
		60020	T-Handle	2	May be used on Quick Angle Handle to angle blade 45° up or down.
		50000G	Instrument Case (23" x 11" x 3 1/2")	1	Holds instruments.





## Urology Blade Kit



### UROLOGY BLADE KIT #SL82064

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
		SL46116	Malleable 13mm x 152mm (½" x 6")	1	Soft tissue retraction.
		SL46122	Malleable 51mm x 203mm (2" x 8")	1	
		SL46123	Malleable 38mm x 203mm (1 ½" x 8")	1	
		SL46124	Malleable 64mm x 254mm (2 ½" x 10")	1	
		SO46129	Kelly 64mm x 76mm (2 ½" x 3")	2	Abdominal or incisional wall retraction.
		SO46136	Fence 102mm x 152mm (4" x 6")	1	Lateral retraction of abdominal wall on large procedures.
		SL46137	Holtgrewe 60mm x 152mm (2 ¾" x 6")	1	Malleable bladder retractor, designed to protect bladder during resection of prostate gland.
		SL46140P	Balfour [with lips] 76mm x 51mm (3" x 2")	2	Costal margin or abdominal wall retraction.
		SL46142P	Balfour [with lips] 51mm x 54mm (2" x 2")	2	
		SL46143	Fenestrated 44mm x 89mm (1 ¾" x 3 ½")	2	Abdominal or incisional wall retraction.
		SO46151	Richardson 51mm x 127mm (2" x 5")	1	Abdominal wall retraction.

UROLOGY BLADE KIT #SL82064 ( Continued )

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
		SL46160	Harrington 64mm x 152mm (2 1/2" x 6")	1	Organ or soft tissue retraction.
		SL46163	Harrington 64mm x 203mm (2 1/2" x 8")	1	
		SL46174	Deaver 51mm x 127mm (2" x 5")	1	Soft tissue retraction.
		SL45753AM	Prostate Retractor 8mm (3/8")	1	Prostate retraction.
		SL46660	Double Angled 25mm x 76mm (1" x 3")	2	Soft tissue retraction.
		SL46661	Double Angled 25mm x 64mm (1" x 2 1/2")	2	
		SL46662	Double Angled 13mm x 76mm (1/2" x 3")	2	
		SL46663	Double Angled 13mm x 70mm (1/2" x 2 3/4")	2	
		SL46215	Murphy Blunt 25mm x 13mm (1" x 1/2")	1	Skin and shallow muscle retraction.
		50000G	Instrument Case (22" x 11" x 3 1/2")	1	Holds instruments.