


Trauma RingTrack® Retractor Set

Abdominal Trauma Set Up



RINGTRACK® FRAME KIT #SL91081

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
Tray 1					
A		47902ACL	RingTrack Elite III Rail Clamp with 1 Cam Joint 16"	1	Attaches to the rail of the OR table and holds the crossbar in position.
B		47437	RingTrack Horizontal Bar 14 1/4"	1	Attaches to the rail clamp and holds the ratcheted ring in place.
C		47412	RingTrack Oval Ring Medium 10 1/2 x 14 1/2"	1	Ratcheted Ring that holds retractor handles, ratchets, and blades around the incision.
D		SL45005G	RingTrack Quick Angle Handle 8"	6	Attaches to S-Lock retractor blades, and inserts into a ratchet to be attached to ring.
E		47440	RingTrack Ratchet	6	Attaches handles to frame.
		60020	T-Handle	2	May be used on Quick Angle Handle to angle blade 45° up or down.
		50000W	Instrument Case (24" x 16" x 3")	1	Holds instruments.

TRAUMA BLADE KIT #SL82052

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
Tray 1					
		SL46112	Malleable 51mm x 165mm (2" x 6 1/2")	2	Soft tissue retraction.
		SL46120B	Malleable 38mm x 152mm (1 1/2" x 6")	1	
		SO46124A	Malleable 76mm x 254mm (3" x 10")	1	
		SO46129	Kelly 64mm x 76mm (2 1/2" x 3")	2	Abdominal or incisional wall retraction.
		SO46140C	Balfour [sm lips] 66mm x 73mm (2 1/2" x 3")	2	Costal margin or abdominal wall retraction.
		SL46143	Fenestrated 44mm x 89mm (1 3/4" x 3 1/2")	2	Abdominal or incisional wall retraction.
		SL46150	Richardson 51mm x 102mm (2" x 4")	1	Abdominal wall retraction.
		SL46151	Richardson 51mm x 127mm (2" x 5")	1	
		SL46152	Richardson 51mm x 152mm (2" x 6")	1	
		SL46160	Harrington 64mm x 152mm (2 1/2" x 6")	1	Organ or soft tissue retraction.
		SO46188	St. Marks [lip] 64mm x 178mm (2 1/2" x 7")	1	Soft tissue retraction.
		50000G	Instrument Case (23" x 11" x 3 1/2")	1	Holds instruments.