


# Bilateral GYN Oncology

















## GYN ONCOLOGY SYSTEM #SL82011

Non S-Lock system also available, may require a longer lead time.

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
<b>Tray 1</b>					
A		43903AC	Elite III Rail Clamp with 2 Cam Joints 18"□	2	Attaches to OR table, holds crossbar, arms, bar.
		41900F	Bilateral Crossbar 7 ¼" x 13 ½" x 7 ¼"	1	Attaches to the rail clamps to provide bilateral stability, and holds the retractor handles in place around the incision.
B		44112SP	Lower Abdominal Bar 14" x 6"	1	Attaches to the rail clamps or crossbar and holds the retractor handles in place around the incision.
		43900BC	24" Crossbar with 2 Cam Joints	1	Attaches to the rail clamp and holds the angled arms or lower abdominal bar in place around the incision. May also support retractor handles.
		50000BL	Instrument Case (26" x 10" x 5")	1	Holds instruments.
<b>Tray 2</b>					
C		41900PH	34" Bilateral Crossbar Hinged (11 ½" x 8 ½" x 11 ½")	1	Attaches to the rail clamps to provide bilateral stability, and holds the retractor handles in place around the incision.
		43112C	Cam III Serrated Slide-On Joint	1	Connects angled arms together to allow for expansion of frame.
		44118N	18" Angled Arm (8" x 10") @ 45°	2	Attaches to the crossbar and holds the retractor handles in place around the incision.
		44124W	24" Angled Arm (11" x 13") @ 45°	2	
D		SL42126WGP	Cam II Clip-on Quick Angle 8"	6	Attaches retractor blades to frame and angles retractor blades as needed, 45° up or down.
E		SL45006CGP	Micro-Adjustable II Clip-on Quick Angle 10"	2	Attaches blades to frame and angles blades as needed, 45° up or down. Provides precise micro-adjustments.
		SL45006CGRP	Micro-Adjustable II Clip-on Quick Angle 15"	2	
		60020	T-Handle	2	Use on QA Handle to angle blade 45° up/down.
		50000G	Instrument Case (23" x 11" x 3 ½")	1	Holds instruments.

**GYN ONCOLOGY SYSTEM #SL82011 (Continued)**

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
<b>Tray 3</b>					
		SL46127	Kelly 38mm x 51mm (1 1/2" x 2")	1	Abdominal or incisional wall retraction.
		SL46128	Kelly 51mm x 64mm (2" x 2 1/2")	1	
		SO46129	Kelly 64mm x 76mm (2 1/2" x 3")	1	
		SO46130	Kelly 76mm x 89mm (3 1/2" x 3")	2	
		SL46115	Malleable 13mm x 64mm (1/2" x 2 1/2")	1	Soft tissue retraction.
		SL46116	Malleable 13mm x 152mm (1/2" x 6")	1	
		SL46117	Malleable 13mm x 203mm (1/2" x 8")	2	
		SL46118	Malleable 25mm x 203mm (1" x 8")	1	
		SL46118A	Malleable 25mm x 254mm (1" x 10")	1	
		SL46118P	Soft Malleable 25mm x 102mm (1" x 4")	1	
		SL46122	Malleable 51mm x 203mm (2" x 8")	1	
		SL46124	Malleable 64mm x 254mm (2 1/2" x 10")	2	
		SO46124A	Malleable 76mm x 254mm (3" x 10")	1	
		SO46124B	Malleable 102mm x 254mm (4" x 10")	1	
		SO46135	Fence 102mm x 127mm (4" x 5")	1	Lateral retraction of abdominal wall on large procedures.
		SO46140C	Balfour [sm lips] 66mm x 73mm (2 1/2" x 3")	1	Costal margin or abdominal wall retraction.
		SO46141C	Balfour [sm lips] 78mm x 73mm (3" x 3")	2	
		SO46141D	Balfour [with lips] 83mm x 123mm (3 1/4" x 5")	2	
		SL46142	Balfour [with lips] 32mm x 32mm (1 1/4" x 1 1/4")	2	
		SL46143	Fenestrated 44mm x 89mm (1 3/4" x 3 1/2")	2	Abdominal or incisional wall retraction.
		SL46150	Richardson 51mm x 102mm (2" x 4")	1	Abdominal wall retraction.
		SO46151	Richardson 51mm x 127mm (2" x 5")	1	
		SO46153	Richardson 57mm x 178mm (2" x 7")	1	
		SL46160	Harrington 64mm x 152mm (2 1/2" x 6")	1	Organ or soft tissue retraction.
		SO46170B	Non-Slip Balfour R 77mm x 66mm (3" x 2 1/2")	1	Costal margin retraction. Aligns with anatomy, reducing slippage.
		SO46171B	Non-Slip Balfour L 77mm x 66mm (3" x 2 1/2")	1	
		SL46174	Deaver 51mm x 127mm (2" x 5")	1	Soft tissue retraction.
		SL46180	Heany 25mm x 102mm (1" x 4")	1	Abdominal or incisional wall retraction.
		SO46188	St. Mark's 64mm x 178mm (2 1/2" x 7")	1	Soft tissue retraction.
		SL46191	Malleable Finger 152mm (6")	2	Gentle organ retraction.
		SO46241A	Wide Balfour 83mm x 57mm (3 1/4" x 2 1/4")	1	Costal margin or abdominal wall retraction.
		SO46241B	Wide Balfour 100mm x 61mm (4" x 2 3/8")	1	
		SO46241C	Wide Balfour 114mm x 63mm (4 1/2" x 2 1/2")	2	
		SL46660	Double Angled 25mm x 76mm (1" x 3")	2	Soft tissue retraction.
		SL46663	Double Angled 13mm x 70mm (1/2" x 2 3/4")	2	
		50000G	Instrument Case (23" x 11" x 3 1/2")	1	Holds instruments.