

## Liver / Oncology System Trays

### Keep it Organized

Thompson's Liver / Oncology System ( #SL82007 ) \* contains three instrument trays. This guide outlines which components to place in each tray to remain within the 25 lb limit.



### LARGE INSTRUMENT CASE 50000BL

50000BL (26" x 10" x 5") can contain both Elites, bilateral crossbar, and Lower Abdominal Bar to equal 22.10 lbs total.

QTY	ITEM DESCRIPTION	PART #
2	Elite III Rail Clamp with 2 Joints 18"	43903AC
1	34" Bilateral Crossbar Hinged (11 1/2" x 8" x 11 1/2")	41900PH
1	Bilateral Crossbar 7 1/4" x 13 1/2" x 7 1/4"	41900F
1	Lower Abdominal Bar 20" (14" x 6")	44112SP

### STANDARD INSTRUMENT CASE 50000G ( 1 OF 2 )

One 50000G (23" x 11" x 3 1/2") can hold the rest of the frame components including handles, joints, and T-Handles to equal 20.48 lbs total. ( Pictured above right. )

QTY	ITEM DESCRIPTION	PART #
1	Angled Arm 18" (8" x 10")	44118N
2	Angled Arm 24" (11" x 13")	44124W
1	Angled Arm 24" (8" x 16")	44124N
2	Cam III Serrated Slide-on Joint	43112C
6	Cam III Clip-on Quick Angle 8"	SL42126WGP
2	Micro-Adjustable II Clip-On Quick Angle 10"	SL45006CGP
2	Micro-Adjustable II Clip-On Quick Angle 15"	SL45006CGRP
2	T-Handle	60020

### STANDARD INSTRUMENT CASE 50000G ( 2 OF 2 )

The second 50000G (23" x 11" x 3 1/2") holds the 37 blades in the system, totaling 13.30 lbs.

QTY	ITEM DESCRIPTION	PART #
2	Malleable 13mm x 64mm (1/2" x 2 1/2")	SL46115
2	Malleable 13mm x 152mm (1/2" x 6")	SL46116
2	Malleable 13mm x 203mm (1/2" x 8")	SL46117
1	Malleable 25mm x 203mm (1" x 8")	SL46118
1	Malleable 25mm x 102mm (1" x 4")	SL46118P
1	Malleable 51mm x 203mm (2" x 8")	SL46122
2	Malleable 64mm x 254mm (2 1/2" x 10")	SL46124
1	Malleable 76mm x 254mm (3" x 10")	SO46124A
1	Malleable 102mm x 254mm (4" x 10")	SO46124B
1	Kelly 64mm x 76mm (2 1/2" x 3")	SO46129
2	Kelly 76mm x 89mm (3" x 3 1/2")	SO46130
1	Balfour [small lips] 70mm x 73mm (2 3/4" x 3")	SO46140C
2	Balfour [with lips] 83mm x 73mm (3 1/4" x 3")	SO46141
2	Balfour [small lips] 83mm x 73mm (3 1/4" x 3")	SO46141C
2	Balfour [with lips] 83mm x 123mm (3 1/4" x 5")	SO46141D
4	Balfour [with lips] 32mm x 32mm (1 1/4" x 1 1/4")	SL46142
1	Richardson 51mm x 178mm (2" x 7")	SL46153
1	Harrington 64mm x 152mm (2 1/2" x 6")	SL46160
1	St. Marks [with lip] 64mm x 178mm (2 1/2" x 7")	SO46188
2	Malleable Finger 152mm (6")	SL46191
1	Wide Balfour 83mm x 57mm (3 3/4" x 2 1/4")	SO46241A
1	Wide Balfour 100mm x 61mm (4" x 2 3/8")	SO46241B
1	Wide Balfour 114mm x 63mm (4 1/2" x 2 1/2")	SO46241C
1	Non-Slip Balfour R 77mm x 94mm (3" x 3 3/4")	SO46169C
1	Non-Slip Balfour L 77mm x 93mm (3" x 3 3/4")	SO46168C

\* Guidelines also apply to Radiolucent Liver / Oncology System ( #SL82008 ).