










# OneFrame™ Frame Kit












## ONEFRAME™ FRAME KIT #SL82063

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
<b>Tray 1</b>					
A		43902AC	Elite III Rail Clamp with 1 Cam Joint 18"	1	Attaches to the rail of the OR table and holds the crossbar in position.
B		47500	OneFrame®	1	Attaches to the rail clamp and holds the retractor handles in place around the incision.
		44013H	Extension Arm 18" with Spring Hinge	1	Attaches to distal end of Angled Armsto provide extra reach, closure, or caudal retraction.
		44118B:	Angled Closure Arm 18" (9" x 9")	2	Attaches to distal ends of frame arms to create enclosed frame.
C		SL42126WGP	Cam II Clip-On Quick Angle 8"	5	Attaches retractor blades to the OneFrame and angles retractor blades as needed, 45° up or down.
		SL45006CGP	Micro-Adjustable II Clip-On Quick Angle 10"	1	Attaches blades to frame and angles blades as needed, 45° up or down. Provides precise micro-adjustments.
		60020	T-Handle	2	May be used on Quick Angle Handle to angle blade 45° up or down.
		50000G	Instrument Case (23" x 11" x 3 1/2")	1	Holds instruments.

**ANTERIOR LUMBAR BLADE KIT #SL82065**

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
<b>Tray 1</b>					
		SL46192ET	Radio. Reverse Lip Malleable 25mm x 200mm	2	Vessel or soft tissue retraction.
		SL46193ET	Radio. Reverse Lip Malleable 32mm x 200mm	2	
		SL46260RET	Radiolucent Reverse Lip 25mm x 100mm	2	Reverse Lip engages at the lateral aspect of the disc and vertebral body to expose the anterior lumbar spine.
		SL46261RET	Radiolucent Reverse Lip 25mm x 110mm	2	
		SL46262RET	Radiolucent Reverse Lip 25mm x 120mm	2	
		SL46263RET	Radiolucent Reverse Lip 25mm x 130mm	2	
		SL46265RET	Radiolucent Reverse Lip 25mm x 150mm	2	
		SL46267RET	Radiolucent Reverse Lip 25mm x 170mm	2	
		SL46269RET	Radiolucent Reverse Lip 25mm x 190mm	2	
		SL46271RET	Radiolucent Reverse Lip 25mm x 210mm	2	
		SL46273RET	Radiolucent Reverse Lip 25mm x 230mm	2	
		SL46275RET	Radiolucent Reverse Lip 25mm x 250mm	2	
		SL46280RET	Radiolucent Reverse Lip 32mm x 100mm	2	
		SL46282RET	Radiolucent Reverse Lip 32mm x 120mm	2	
		SL46284RET	Radiolucent Reverse Lip 32mm x 140mm	2	
		SL46286RET	Radiolucent Reverse Lip 32mm x 160mm	2	
	SL46288RET	Radiolucent Reverse Lip 32mm x 180mm	2		
	SL46290RET	Radiolucent Reverse Lip 32mm x 200mm	2	Reverse Lip engages at the lateral aspect of the disc and vertebral body to expose the anterior lumbar spine.	
	SL46292RET	Radiolucent Reverse Lip 32mm x 220mm	2		
	SL46293RET	Radiolucent Reverse Lip 32mm x 230mm	2		
	SL46294RET	Radiolucent Reverse Lip 32mm x 240mm	2		
		SL46295RET	Radiolucent Reverse Lip 32mm x 250mm	2	
		SL46260TET	Radio. Reverse Lip Tapered 25mm x 100mm	2	Tapered design is ideal for working around osteophytes.
		SL46265TET	Radio. Reverse Lip Tapered 25mm x 150mm	2	
		SL46270TET	Radio. Reverse Lip Tapered 25mm x 200mm	2	
		SL46119BET	Radiolucent Malleable Concave Round Neck Renal 25mm x 140mm	2	Renal vein or vessel retraction.
		SL46119CET	Radiolucent Malleable Concave Round Neck Renal 25mm x 190mm	2	
	( see image below )	SL46564	Radio. Conc. Short Neck Renal 25mm x 140mm	4	Rigid, low profile vessel retraction.
		SL46569	Radio. Conc. Short Neck Renal 25mm x 190mm	4	
		SL46575	Radio. Conc. Short Neck Renal 25mm x 250mm	4	
		51234	Suction for Anterior Lumbar Surgery	1	Suction for disc space, as needed.
		51236	Anterior Lumbar Depth Gauge 335mm	1	Measures depth of incision for blade selection.
		50000ERL	Instrument Case - Exp. Reverse Lip ASA	1	Organizes and holds instruments.
<b>Tray 2</b>					
		SL46558	Radio. Conc. Short Neck Renal 25mm x 80mm	4	Rigid, low profile vessel retraction.
		SL46560	Radio. Conc. Short Neck Renal 25mm x 100mm	4	
		SL46562	Radio. Conc. Short Neck Renal 25mm x 120mm	4	
		SL46564	Radio. Conc. Short Neck Renal 25mm x 140mm	4	
		SL46566	Radio. Conc. Short Neck Renal 25mm x 160mm	4	
		SL46568	Radio. Conc. Short Neck Renal 25mm x 180mm	4	
		SL46570	Radio. Conc. Short Neck Renal 25mm x 200mm	4	
		50000L	Instrument Tray - Lid	1	Tray Lid
		50000RV	Instrument Tray - Renal Vein Blades	1	Drops into Tray Base. Organizes blades.
		50000U1	Instrument Tray - Base	1	Tray Base.
<b>Tray 3</b>					
	( see image above )	SL46119BET	Radio. Malleable Conc. Round Neck Renal 25mm x 140mm	2	Renal vein or vessel retraction.
		SL46119CET	Radio. Malleable Conc. Round Neck Renal 25mm x 190mm	2	
		SL46119RET	Radiolucent Long Neck Renal 25mm x 140mm	2	Renal vein or vessel retraction.
		SL46120RET	Radiolucent Long Neck Renal 25mm x 190mm	2	
		SL46121RET	Radiolucent Long Neck Renal 25mm x 255mm	2	











## ANTERIOR LUMBAR BLADE KIT #SL82065 (Continued)

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
		SL46122ET	Radio. Malleable 51mm x 203mm (2" x 8")	2	Soft tissue retraction.
		SL46140ET	Radio. Balfour 65mm x 72mm (2 1/2" x 2 3/4")	2	Abdominal wall retraction.
		SL46141ET	Radio. Balfour 83mm x 70mm (3 1/4" x 2 3/4")	1	
		SL46145ET	Radio. Splanchnic 51mm x 115mm (2" x 4 1/2")	1	Soft tissue retraction.
		SL46160ET	Radio. Harrington 64mm x 152mm (2 1/2" x 6")	1	Organ or soft tissue retraction.
		SL46176ET	Radiolucent Deaver 76mm x 127mm (3" x 5")	1	Soft tissue retraction.
	( see image above )	51234	Suction for Anterior Lumbar Surgery	1	Suction for disc space, as needed.
		SL42126WDP	Cam II Clip-on Dual Angle 8"	2	Attaches blades to frame and angles blades with tactile positioning and precise angling adjustments.
		SL42128G	Quick Angle Hand Held Adapter	1	Attach to a Thompson retractor blade to use as a hand held retractor.
		60021	Hex Tool 3/16" with T-Handle	2	May be used on Quick Angle and Dual Angle Handles to angle blade up or down.
		50000ALTI	Case - Anterior Lumbar Insert	1	Organizes and Holds instruments
		50000G	Instrument Case (23" x 11" x 3 1/2")	1	Organizes and holds instruments.

## GENERAL / VASCULAR BLADE KIT #SL82066

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
		SL46119B	Concave Round Neck Renal 25mm x 140mm (1" x 5 1/2")	2	Renal vein or vessel retraction.
		SL46119C	Concave Round Neck Renal 25mm x 190mm (1" x 7 1/2")	2	
		SL46122V	Malleable Slotted 51mm x 203mm (2" x 8")	1	Soft tissue retraction.
		SO46124AV	Malleable Slotted 76mm x 203mm (3" x 8")	1	
		SL46118	Malleable 25mm x 203mm (1" x 8")	1	Soft tissue retraction.
		SL46122	Malleable 51mm x 203mm (2" x 8")	2	
		SL46124	Malleable 64mm x 254mm (2 1/2" x 10")	1	
		SO46124A	Malleable 76mm x 254mm (3" x 10")	1	
		SL46125	Richardson 19mm x 25mm (3/4" x 1")	1	Muscle retraction for exposure with small or mini-incision.
		SO46129	Kelly 64mm x 76mm (2 1/2" x 3")	2	Abdominal or incisional wall retraction.
		SO46130	Kelly 76mm x 89mm (3" x 3 1/2")	1	



**GENERAL / VASCULAR BLADE KIT #SL82066 (Continued)**

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
		SO46134	Fence 102mm x 102mm (4" x 4")	1	Lateral retraction of abdominal wall on large procedures.
		SO46135	Fence 102mm x 127mm (4" x 5")	1	
		SO46136	Fence 102mm x 152mm (4" x 6")	1	
		SO46140	Balfour [lips] 70mm x 73mm (2 3/4" x 3")	1	Costal margin or abdominal wall retraction.
		SO46140C	Balfour [sm lips] 66mm x 73mm (2 1/2" x 3")	2	
		SO46141	Balfour [lips] 83mm x 73mm (3 1/4" x 3")	1	
		SO46141C	Balfour [sm lips] 78mm x 73mm (3" x 3")	1	
		SL46143	Fenestrated 44mm x 89mm (1 3/4" x 3 1/2")	2	Abdominal or incisional wall retraction.
		SL46144	Fenestrated 48mm x 127mm (1 7/8" x 5")	2	
		SL46145	Splanchnich 51mm x 114mm (2" x 4 1/2")	1	Soft tissue retraction.
		SL46146	Splanchnich 51mm x 152mm (2" x 6")	1	
		SL46151	Richardson 51mm x 127mm (2" x 5")	2	Abdominal wall retraction.
		SL46160	Harrington 64mm x 152mm (2 1/2" x 6")	1	Organ or soft tissue retraction.
		SL46174	Deaver 51mm x 127mm (2" x 5")	1	Soft tissue retraction.
		SO46176	Deaver 76mm x 127mm (3" x 5")	1	
		SO46188	St. Marks [lip] 64mm x 178mm (2 1/2" x 7")	1	Soft tissue retraction.
		SL46191	Malleable Finger 152mm (6")	1	Gentle organ retraction.
		SO46241A	Wide Balfour 83mm x 57mm (3 1/4" x 2 1/4")	1	Costal margin or abdominal wall retraction.
		SO46241B	Wide Balfour 100mm x 61mm (4" x 2 3/8")	1	
		SO46241C	Wide Balfour 114mm x 63mm (4 1/2" x 2 1/2")	1	
		50000G	Instrument Case (23" x 11" x 3 1/2")	1	Holds instruments.









## UROLOGY BLADE KIT #SL82064

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
		SL46116	Malleable 13mm x 152mm (1/2" x 6")	1	Soft tissue retraction.
		SL46122	Malleable 51mm x 203mm (2" x 8")	1	
		SL46123	Malleable 38mm x 203mm (1 1/2" x 8")	1	
		SL46124	Malleable 64mm x 254mm (2 1/2" x 10")	1	
		SO46129	Kelly 64mm x 76mm (2 1/2" x 3")	2	Abdominal or incisional wall retraction.
		SO46136	Fence 102mm x 152mm (4" x 6")	1	Lateral retraction of abdominal wall on large procedures.
		SL46137	Holtgrewe 60mm x 152mm (2 3/8" x 6")	1	Malleable bladder retractor, designed to protect bladder during resection of prostate gland.
		SL46140P	Balfour [with lips] 76mm x 51mm (3" x 2")	2	Costal margin or abdominal wall retraction.
		SL46142P	Balfour [with lips] 51mm x 54mm (2" x 2")	2	
		SL46143	Fenestrated 44mm x 89mm (1 3/4" x 3 1/2")	2	Abdominal or incisional wall retraction.
		SO46151	Richardson 51mm x 127mm (2" x 5")	1	Abdominal wall retraction.
		SL46160	Harrington 64mm x 152mm (2 1/2" x 6")	1	Organ or soft tissue retraction.
		SL46163	Harrington 64mm x 203mm (2 1/2" x 8")	1	
		SL46174	Deaver 51mm x 127mm (2" x 5")	1	Soft tissue retraction.
		SL45753AM	Prostate Retractor 8mm (3/8")	1	Prostate retraction.





**UROLOGY BLADE KIT #SL82064 (Continued)**

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
		SL46660	Double Angled 25mm x 76mm (1" x 3")	2	Soft tissue retraction.
		SL46661	Double Angled 25mm x 64mm (1" x 2 1/2")	2	
		SL46662	Double Angled 13mm x 76mm (1/2" x 3")	2	
		SL46663	Double Angled 13mm x 70mm (1/2" x 2 3/4")	2	
		SL46215	Murphy Blunt 25mm x 13mm (1" x 1/2")	1	Skin and shallow muscle retraction.
		50000G	Instrument Case (22" x 11" x 3 1/2")	1	Holds instruments.

**KIDNEY BLADE KIT #SL82071**

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
<b>Tray 1</b>					
		SL46119B	Concave Round Neck Renal 25mm x 140mm (1" x 5 1/2")	1	Renal vein or vessel retraction.
		SL46119C	Concave Round Neck Renal 25mm x 190mm (1" x 7 1/2")	1	
		SL46118	Malleable 25mm x 203mm (1" x 8")	1	Soft tissue retraction.
		SL46118A	Malleable 25mm x 254mm (1" x 10")	1	
		SL46118B	Malleable 38mm x 254mm (1 1/2" x 10")	1	
		SL46122	Malleable 51mm x 203mm (2" x 8")	1	
		SL46123	Malleable 38mm x 203mm (1 1/2" x 8")	1	
		SL46124	Malleable 64mm x 254mm (2 1/2" x 10")	1	
		SO46124A	Malleable 76mm x 254mm (3" x 10")	1	
		SL46127	Kelly 38mm x 51mm (1 1/2" x 2")	1	Abdominal or incisional wall retraction.
		SL46128	Kelly 51mm x 64mm (1 1/2" x 2 1/2")	1	
		SO46129	Kelly 64mm x 76mm (2 1/2" x 3")	1	
		SL46140P	Balfour 76mm x 51mm (3" x 2")	1	Costal margin or abdominal wall retraction.
		SL46142P	Balfour 51mm x 54mm (2" x 2")	1	
		SL46125	Richardson 19mm x 25mm (3/4" x 1")	1	Muscle retraction for exposure with small or mini-incision.
		SL46125A	Richardson 19mm x 51mm (3/4" x 2")	1	
		SL46126	Richardson 25mm x 32mm (1" x 1 1/4")	1	
		SL46150	Richardson 51mm x 102mm (2" x 4")	1	Abdominal wall retraction.
		SO46151	Richardson 51mm x 127mm (2" x 5")	1	
		SL46154	Richardson 51mm x 76mm (2" x 3")	1	
		SL46155	Mini Open Richardson 51mm x 76mm (2" x 3")	1	Muscle retraction for exposures with small or mini-incision.
		SL46156	Mini Open Richardson 51mm x 102mm (2" x 4")	1	
		SL46157	Mini Open Richardson 51mm x 127mm (2" x 5")	1	
		SL46158	Mini Open Richardson 51mm x 152mm (2" x 6")	1	
		SL46160	Harrington 64mm x 152mm (2 1/2" x 6")	1	Organ or soft tissue retraction.
		SL46162P	Harrington 38mm x 127mm (1 1/2" x 5")	1	
		SL46163	Harrington 64mm x 203mm (2 1/2" x 8")	1	
		SL46165	Harrington 64mm x 102mm (2 1/2" x 4")	1	

**KIDNEY BLADE KIT #SL82071 (Continued)**

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
		SL46220	Mini-Open Malleable 25mm x 75mm (1" x 3")	1	Soft tissue retraction for exposures with small or mini-incision.
		SL46221	Mini-Open Malleable 25mm x 100mm (1" x 4")	1	
		SL46222	Mini-Open Malleable 25mm x 125mm (1" x 5")	1	
		SL46224	Mini-Open Double Angled 51mm x 127mm (2" x 5")	1	Muscle retraction for exposures with small or mini-incision.
		SL46225	Mini-Open Double Angled 51mm x 152mm (2" x 6")	1	
		SL46226	Mini-Open Double Angled 51mm x 165mm (2" x 6 1/2")	1	
		SL46227	Mini Open Deep Double Angled 51mm x 127mm (2" x 5")	1	Muscle retraction for exposures with small or mini-incision.
		SL46660	Double Angled 25mm x 76mm (1" x 3")	1	Soft tissue retraction.
		50000G	Instrument Case (23" x 11" x 3 1/2")	1	Holds instruments.