







Open Posterior Lumbar

















POSTERIOR LUMBAR SYSTEM #SL82020

Non S-Lock system also available, may require a longer lead time.

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
Tray 1					
A		43902ACL	Elite III Rail Clamp with 1 Cam Joint 16"	1	Attaches to the rail of the OR table and holds the crossbar in position.
B		43900BC	Crossbar w/2 Cam Joints 24" (11" x 13")	1	Attaches to the rail clamp and holds the angled arms in place around the incision. May also support retractor handles.
C		44124N	24" Angled Arm (8" x 16") @ 45°	2	Attaches to the crossbar and holds the retractor handles in place around the incision.
D		42140C	Flexible Blade Holder with Cam Joint	1	Attaches Nerve Root Retractors to the frame as needed.
E		SL45011CGP	Mini Micro-Adjustable II Clip-on Quick Angle 9"	4	Attaches blades to frame and angles blades as needed, 45° up or down. Provides precise micro-adjustments. Lower profile.

POSTERIOR LUMBAR SYSTEM #SL82020 (Continued)

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
		60020	T-Handle	2	May be used on Quick Angle Handle to angle blade 45° up or down.
		50000G	Instrument Case (23" x 11" x 3 1/2")	1	Holds instruments for sterilization, storage, and transportation.
Tray 2					
		SO45670ATTC	Radiolucent Concave Posterior Lumbar 23mm x 80mm (1" x 3 1/8")	2	Lateral muscle retraction for posterior lumbar procedures.
		SO45671ATTC	Radiolucent Concave Posterior Lumbar 45mm x 80mm (1 3/4" x 3 1/8")	2	
		SO45672ATTC	Radiolucent Concave Posterior Lumbar 60mm x 80mm (2 3/8" x 3 1/8")	2	
		SO45673ATTC	Radiolucent Concave Posterior Lumbar 23mm x 100mm (1" x 4")	2	
		SO45674ATTC	Radiolucent Concave Posterior Lumbar 45mm x 100mm (1 3/4" x 4")	2	
		SO45675ATTC	Radiolucent Concave Posterior Lumbar 60mm x 100mm (2 3/8" x 4")	2	
		SO45676ATTC	Radiolucent Concave Posterior Lumbar 23mm x 120mm (1" x 4 3/4")	2	
		SO45677ATTC	Radiolucent Concave Posterior Lumbar 45mm x 120mm (1 3/4" x 4 3/4")	2	
		SO45678ATTC	Radiolucent Concave Posterior Lumbar 60mm x 120mm (2 3/8" x 4 3/4")	2	
		45752A	Nerve Root Retractor 90° 5mm x 123mm (1/4" x 4 7/8")	1	Held in place by the Flexible Blade Holder and provide posterior lumbar nerve root retraction.
		45753A	Nerve Root Retractor 45° 8mm x 145mm (3/8" x 5 3/4")	1	
		SL45760A	Prong Retractor 75mm (3")	1	Retracts off of the spinous process, or bone, to increase access around the retractor for procedures such as drilling
		SL45761A	Prong Retractor 100mm (4")	1	
		50000G	Instrument Case (23" x 11" x 3 1/2")	1	Holds instruments.