




Hip



Knee





Shoulder



TOTAL JOINT FRAME #SL82045

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
Tray 1					
A		43904AC	Elite III Rail Clamp with 1 Joint 22"	1	Attaches to the rail of the OR table and holds the crossbar in position.
B		43990	Crossbar with 2 Cam III Joints 7" x 19 1/2"	1	Attaches to the rail clamp and holds the curved arms in place around the incision. May also support retractor handles.
C		44203	20" Curved Arm with Cleats 1	1	Attaches to the crossbar and holds the retractor handles and tethers in place around the incision.
D		44204	20" Curved Arm with Cleats 2	1	
E		42105	Hip Tethers 24" - 5 Pack (Latex Free)	2	Secures retractor blades to the frame.
F		42126BP	Cam II Clip-On Blade Adapter 8 1/2"	4	
		42126OHP	Orthopedic Cam II Clip-On Hook 8 1/2"	1	Secures retractor blades to the frame.

TOTAL JOINT FRAME #SL82045 (Continued)

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
		41938	Rail Extender 20" with 1 Power Rail Clamp	1	Attaches to OR table rail to provide more attachment options for the rail clamp. Required to use Hip System with Hana® Table or leg positioning systems.
		50000G	Instrument Case (23" x 11" x 3 1/2")	1	Holds instruments for sterilization, storage, and transportation.

HIP BLADE KIT #82047

REF	IMAGE	PART #	ITEM DESCRIPTION	QTY	APPLICATION
		45301	Cobra Hip Blade 11"	2	
		45305	Hohmann Hip Blade 11"	2	Acetabulum exposure.
		45309	Dual Prong Hip Blade 11"	2	
		45300	Femoral Neck Elevator 14"	1	Femur elevation.
		45302	Curved Cobra Hip Blade 12 1/2"	1	
		45306	Curved Hohmann Narrow Hip Blade 12"	1	
		45307	Curved Hohmann Wide Hip Blade 12"	1	Acetabulum exposure.
		45308	Curved Dual Prong Hip Blade 12"	1	
		50000G	Instrument Case (23" x 11" x 3 1/2")	1	Holds instruments for sterilization, storage, and transportation.