



## Perineal Prostatectomy



## Retropubic Prostatectomy













### UROLOGY SYSTEM #SL82012

Non S-Lock system also available, may require a longer lead time.

| REF           | IMAGE | PART #     | ITEM DESCRIPTION                            | QTY | APPLICATION  |
|---------------|-------|------------|---|-----|--|
| <b>Tray 1</b> |       |            |   |     |  |
| A             |       | 43902AC    | Elite III Rail Clamp with 1 Cam Joint 18"   | 1   | Attaches to the rail of the OR table and holds the crossbar in position.   |
| B             |       | 43900BC    | 24" Crossbar with 2 Cam III Joints          | 1   | Attaches to the rail clamp and holds the angled arms in place around the incision. May also support retractor handles. |
| C             |       | 44118N     | 18" Angled Arm (8" x 10") @ 45°             | 2   | Attaches to the crossbar and holds the retractor handles in place around the incision.                                 |
|               |       | 44124N     | 24" Angled Arm (8" x 16") @ 45°             | 2   |  |
| D             |       | SL42126WGP | Cam II Clip-on Quick Angle 8"               | 3   | Attaches retractor blades to the crossbar or angled arm and angles retractor blades as needed, 45° up or down.         |
|               |       | SL45006CGP | Micro-Adjustable II Clip-On Quick Angle 10" | 2   |  |
|               |       | 60020      | T-Handle                                    | 2   | Use on QA Handle to angle blade 45° up/down.   |
|               |       | 50000G     | Instrument Case (23" x 11" x 3 1/2")        | 1   | Holds instruments.   |
| <b>Tray 2</b> |       |            |   |     |  |
|               |       | SL46116    | Malleable 13mm x 152mm (1/2" x 6")          | 1   | Soft tissue retraction.  |
|               |       | SL46122    | Malleable 51mm x 203mm (2" x 8")            | 1   |  |
|               |       | SL46123    | Malleable 38mm x 203mm (1 1/2" x 8")        | 1   |  |
|               |       | SL46124    | Malleable 64mm x 254mm (2 1/2" x 10")       | 1   |  |

UROLOGY SYSTEM #SL82012 (Continued)

| REF | IMAGE   | PART #    | ITEM DESCRIPTION                          | QTY | APPLICATION  |
|-----|---|-----------|---|-----|--|
|     |    | SO46129   | Kelly 64mm x 76mm (2 1/2" x 3")           | 2   | Abdominal or incisional wall retraction.   |
|     |    | SO46136   | Fence 102mm x 152mm (4" x 6")             | 1   | Lateral retraction of abdominal wall on large procedures.                                    |
|     |    | SL46137   | Holtgrewe 60mm x 152mm (2 3/8" x 6")      | 1   | Malleable bladder retractor, designed to protect bladder during resection of prostate gland. |
|     |    | SL46140P  | Balfour [with lips] 76mm x 51mm (3" x 2") | 2   | Costal margin or abdominal wall retraction.  |
|     |   | SL46142P  | Balfour [with lips] 51mm x 54mm (2" x 2") | 2   |  |
|     |    | SL46143   | Fenestrated 44mm x 89mm (1 3/4" x 3 1/2") | 2   | Abdominal or incisional wall retraction.   |
|     |   | SO46151   | Richardson 51mm x 127mm (2" x 5")         | 1   | Abdominal wall retraction.   |
|     |  | SL46160   | Harrington 64mm x 152mm (2 1/2" x 6")     | 1   | Organ or soft tissue retraction.   |
|     |   | SL46163   | Harrington 64mm x 203mm (2 1/2" x 8")     | 1   |  |
|     |  | SL46174   | Deaver 51mm x 127mm (2" x 5")             | 1   | Soft tissue retraction.  |
|     |  | SL45753AM | Prostate Retractor 8mm (3/8")             | 1   | Prostate retraction.   |
|     |  | SL46660   | Double Angled 25mm x 76mm (1" x 3")       | 2   | Soft tissue retraction.  |
|     |   | SL46661   | Double Angled 25mm x 64mm (1" x 2 1/2")   | 2   |  |
|     |   | SL46662   | Double Angled 13mm x 76mm (1/2" x 3")     | 2   |  |
|     |   | SL46663   | Double Angled 13mm x 70mm (1/2" x 2 3/4") | 2   |  |
|     |  | SL46215   | Murphy Blunt 25mm x 13mm (1" x 1/2")      | 1   | Skin and shallow muscle retraction.  |
|     |   | 50000G    | Instrument Case (22" x 11" x 3 1/2")      | 1   | Holds instruments.   |